

5ml & 10ml Pipette Tips

I. APPLICATIONS

1. For use with pipettors for the measure and transfer of liquids/chemicals
2. Fit Thermo 5ml/10ml Pipette

II. FEATURES

1. Made of medical-grade PP material, allowing for easy observation of the liquid status in the pipette tip.
2. Chemically resistant to alcohols and mild organic solvents.
3. 100% inert.
4. Certified RNase & DNase Free, non-pyrogen.
5. Perfectly adapt with Thermo and Brand pipettes.

III. FILTER TIPS DESIGN

1. 5ml and 10ml tip sizes.
2. High sealing performance to prevent air leakage or back-suction.
3. Clearly marked volume scale for quick observation of the pipetting capacity, allowing for calibration and effectively avoiding errors.
4. Design of smooth inner surface effectively prevents any liquid leakage or residue, ensuring smooth and efficient pipetting.

IV. DESCRIPTION

Fit Thermo 5ml/10ml Pipette

Cat. #	Size	With Filter	Sterile	Packaging Type	Packaging	Qty/Case
21-5000TS	5ml	No	Yes	Rack Pack	40/Rack, 40 Racks/Case	1600
20-5000TS	5ml	No	No	Bulk Pack	100/Bag, 20 Bags/Case	2000
21-1100TS	10ml	No	Yes	Rack Pack	24/Rack, 40 Racks/Case	960
20-1100TS	10ml	No	No	Bulk Pack	100/Bag, 10 Bags/Case	1000
23-5000TS	5ml	Yes	Yes	Rack Pack	40/Rack, 40 Racks/Case	1600
22-5000TS	5ml	Yes	No	Bulk Pack	100/Bag, 20 Bags/Case	2000
23-1100TS	10ml	Yes	Yes	Rack Pack	24/Rack, 40 Racks/Case	960
22-1100TS	10ml	Yes	No	Bulk Pack	100/Bag, 10 Bags/Case	1000

V. STERILITY

1. Yes, for 21- & 23-

VI. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2015 standard.

VII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

VIII. Photos



IX. PRODUCT PARAMETER

Dimensions=mm

