

Petri Dishes

I. APPLICATIONS

1. Microbiology
2. Cell growth

II. FEATURES

1. Materials: polystyrene
2. EO sterilized
3. Disposable
4. Glass-like clarity

III. PETRI DISH DESIGN

1. Available in 150 x 15mm, 90 x 15mm or 60 x 15mm
2. Triple-vented lid
3. Stackable

IV. PACKAGING

1. Bulk, 500 pieces/case (90 x 15mm or 60 x 15mm)
2. Bulk, 200 pieces/case (150 x 15mm)

V. DESCRIPTION

Item	Description	Packaging	Sterile
66-1515	150 x 15mm petri dish	Bulk, 200 pieces/case	Yes
66-1501	90 x 15mm petri dish	Bulk, 500 pieces/case	Yes
66-1560	60 x 15mm petri dish	Bulk, 500 pieces/case	Yes

VI. STERILITY

1. Yes.
2. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2015 standard.





VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H
66-1501	Petri Dish	O.D. 89.3	O.D. 86.5	7.9	6.7	14.9	13.4	16.1	1.2
66-1560	Petri Dish	O.D. 60	O.D. 57	7	/	13	/	15.5	/
66-1515	Petri Dish	O.D. 150	O.D. 144	8.7	6.6	14.3	12.2	16.6	2.1



66-1501

66-1560

